Ethiestic Plus

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Ethiestic Plus

Silicone brushes for shaping composites

Ethiestic Plus were developed to promote easiness and precision in shaping composites, both direct and indirect restorations. The kit contains a straight instrument, which facilitates the access to the vestibular anterior/posterior teeth, as well as for laboratory work, and an angled instrument, which allows an easier and more comfortable access to the lingual/palatal region of the anterior/posterior teeth and to the occlusal region of posterior teeth. The brushes are interchangeable, made of silicone rubber, anti-adherent and available in four different shapes and two different hardness (hard and soft tips), which allows to shape the composite to obtain the desired tooth anatomy. The possibility of rotating the brushes on the end of the instrument provides the comfort during its use. Using Ethiestic Plus avoids the typical lines and marks caused by hard instruments and, consequently, the differences of shine and color that they may cause.

Use: Select the instrument and brush (tip), which best fits to the clinical case. Place the brush on the end of the instrument. Ensure that the brush is completely fitted in the instrument. Rotate the brush to reach the best position. Shape the composite with light touches, starting with the hard brushes (blue) and finishing with the soft brushes (white).

Warnings, precautions and restrictions: Professional use. If there is proven allergy to any of its components, the product must not be used. To avoid an accidental swallowing of the product, it is necessary the application of a total isolation technique of the operating field. Product not recommended for: condensable composites, flow composites, amalgams, glass ionomer cements and self-curing acrylic resins. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Cleaning: Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amylase, Protease, Lipase and Carboxylase) and with pure pH range between 6.5 - 7.5. Prepare the dilution solution according to the manufacturer's instructions. The water should be potable. For best effectiveness, it is recommended that the temperature be between 40°C and 55°C. Submerge the product and leave in contact for 5 minutes with the solution. Remove parts and rinse under running water for at least 10 seconds. Dry the parts with absorbent paper or disposable wipes. **Sterilization:** The instrument should be autoclaved (126°C, 16 min.), after each patient. Sterilization performed with any method other than autoclave can damage the product. Sterilization methods may cause a slight color changing in the brushes, which does not affect its functionality. **Inspection:** Observe damage / breakage of the product that would prevent proper operation and do not use in this case. **Reprocessing Limits:** This is a reusable device. Repetitive processing has minimal effect due to its 200 cycles for the instrument, and 100 cycles for the brushes. The end of the life, identified by the professional, is determined due to damage resulting from use.

Storage: To be stored in a cool and fresh place, protected from sunlight.

Dimensions: Available in 4 different shapes: Double chamfered angled set, single chamfered angled set, double angled chamfered angled set, flame.

Composition: Instruments: Polysulfone. Soft brushes: Silicone and titanium dioxide. Hard brushes: Silicone, titanium dioxide, aluminum oxide and atoxic colorants.

Disposal: Dispose in accordance with applicable local requirements.

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Manufacturer: TDV Dental Ltd.
EC REP/OBELIS S.A.