Unimatrix R

Hard ◯ Soft

- PT Sistema de matrizes seccionais pré-formadas para restaurações de Classe II
- ES Sistema de matrices seccionales pre-formadas para restauraciones Clase II
- EN Preformed sectional matrices system for Class II restorations
- FR Système de matrices sectionnelles pour des restaurations de Classe II
- DE System für vorgeformte Teilmatrizen für Restaurationen der Klasse II
- IT Sistema di matrice sezionale preformata per i restauri di Classe II
- NL Systeem van voorgevormde sectionele matten voor klasse II restauraties
- CS Systém práce s sekčními matricemi u oprav třídy II

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INSTRUCTIONS FOR USE (EN)

The Unimatrix R System is indicated for Class II restorations and composed by preformed sectional matrices, available in 3 sizes, and reinforced rings. The main advantages of this system are the easiness in the reconstruction of contact points and natural tooth anatomy. The reinforced rings are equipped with triangular silicone tips, and available in two models: Hard (blue) indicated for clinical cases requiring more pressure and dental separation, and Soft (pink) indicated for clinical cases requiring less pressure and dental separation.

Uso: Isolate the operative field with a rubber dam or other suitable isolation technique. Select, among the 3 matrix sizes, the item which is best adapted to the cavity to be restored. Install the matrix in the proximal surface to be restored. When it concerns a MOD cavity, it is advisable to restore one surface at a time. Install the wedge. TDV elastic wedges improve the adaptation of the matrix. In this case, silicone tips are not used and the ring is positioned in the elastic wedge piercings. Install the Soft (pink) or Hard (blue) Unimatrix R ring as per the clinical case. For opening and inserting the ring, TDV Unimatrix pliers or Palmer pinchers can be used. The use of or not of the ring and silicone tips and the means of their installation depend upon personal preferences and the clinical situation. Having concluded the restoration, remove first the fastening clamp, then the wedge and lastly the matrix. To obtain consistent contact points, it is advisable to separate the adjacent ones by means of wedges or fastening ring, from the anesthetic until the installation of the matrix.

Warnings, precautions and restrictions: Professional use. The matrices are of single use. Their reuse can cause cross-contamination. When there is proven allergy to any one of its components, the product must not be used. Excessive expansion with the pliers can lead to deformation of the ring. During the insertion of the ring, maintain the finger pressed on it until the removal of the pliers, ensuring the adjustment and avoiding sliding. Care must be taken during the insertion of the ring, due to its flexibility. Ensure that the product is totally fitted in the pliers. Avoid leaving the ring open in the pinchers for a long period of time. It is advisable to protect the papilla by any method to avoid accidents. It is advisable that the use of the product is always accompanied by absolute isolation of the operative field, as objects can be swallowed or breathed in by accident. The ring must be used with the silicone protectors for an appropriate adjustment of the matrix. The ring and the silicone tips must be sterilized separately. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

RING / SILICONE TIPS - Cleaning: Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amylase, Protease, Lipase and Carbohydrase) and with pure pH range between 6.5 - 7.5. Prepare the solution according to the manufacturer's instructions. The water should be potable. For a better effectiveness, it is recommended that the temperature be between 40°C and 55°C. Submerge the products and leave in contact for 5 minutes with the solution. Remove parts and rinse under running water for at least 10 seconds. Dry the pieces with absorbent paper or disposable wipes. **Sterilization:** The ring and the triangular silicone tips must be sterilized separately in autoclave until 134 °C - 6 min. The sterilization can cause slight change in the coloring of the tips and reinforcement, but does not affect its functionality. Check the product before each use. Note if there is damage (breaks/cracks/rust) that prevents proper operation and does not use in such cases. These conditions are indicative of fragility of the material, whose strength can be reduced during use. **Reprocessing Limits:** This is a reusable device. Repeated processing has minimal effect due to its composition. The end of the life of the ring and triangular silicone tips, identified by the professional, is determined due to damage resulting from use.

Composition: Matrices: stainless steel. Ring: nickel-plated carbon steel. Ring reinforcement: thermoplastic. Tips: silicone

Storage: To be stored in a cool and dry place, protected from sunlight.

Disposal: Dispose in accordance with applicable local requirements

Ref.	Code	
04126H	0003020081	Unimatrix R ring Hard – 1 ring + 4 triangular silicone tips
04124R	0003020085	Unimatrix R – 50 matrices (S/M/L) + 2 rings + 8 triangular silicone tips
04122R	0003020086	Unimatrix R – 25 matrices (S/M/L) + 1 ring + 4 triangular silicone tips
04120R	0003020087	Unimatrix R – 100 matrices (S/M/L) + 2 rings + 8 triangular silicone tips
04126S	0003020102	Unimatrix R ring Soft – 1 ring + 4 triangular silicone tips
04159T	0003020103	Triangular silicone tips for Unimatrix rings – 20 units
04160T	0003020104	Triangular silicone tips for Unimatrix rings – 50 units
-	0003020105	Unimatrix R – 100 Matrices (S/M/L) + 2 rings (Hard/Soft) + 50 triangular silicone tips + 50 elastic wedges
-	0003020106	Unimatrix R – 50 Matrices (S/M/L) + 2 rings (Hard/Soft) + 20 triangular silicone tips + 25 elastic wedges

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Obelis S.A